


Product Datasheet

Rb (phospho-Ser807) antibody (orb559130)

Description	Rb (phospho-Ser807) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Rb around the phosphorylation site of Ser807 (YI(p-S)PLK)
Target	Phospho-Rb (Ser807)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	102 kDa
Uniprot ID	P06400
Dilution Range	IHC-P=1:50-200, WB=1:500-2000, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
colon
carcinoma
...