

Product Datasheet

PLK1 (phospho-Thr210) antibody (orb559133)

Description	PLK1 (phospho-Thr210) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PLK1 around the phosphorylation site of Thr210 (KK(p-T)L)
Target	Phospho-PLK1 (Thr210)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	68 kDa
CAS Number	Serine/threonine-protein kinase PLK1
Uniprot ID	P53350
Dilution Range	Flow-Cyt=1:50, WB=1:500-1000, IHC-P=1:100-500, IHC-F=1:10-50, ICC=1:10-50, IF=1:10-50
Expiration Date	12 months from date of receipt.