

Product Datasheet

PKR (phospho-Thr446) antibody (orb559134)

Description

PKR (phospho-Thr446) antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PKR around the phosphorylation site of Thr446 (KR(p-T)RS)

Target

Phospho-PKR (Thr446)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Recombinant

MW

62 kDa

Uniprot ID


P19525

Dilution Range


WB=1:500-1000, IHC-P=1:50-200

Expiration Date


12 months from date of receipt.



IHC-P
analysis of
mouse
placenta
using P...



IHC-P
analysis of
rat brain
using
Phosph...



IHC-P
analysis of
human
colon
carcinoma
...