



Product Datasheet

AGR2 antibody (orb559155)

Description	AGR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human AGR2
Target	AGR2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	20 kDa
Uniprot ID	O95994
Dilution Range	WB=1:1000-5000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human lung
carcinoma
u...



Western
blot: Lane
1: A549
(Human)
Cell; Pr...