


Product Datasheet

P53 (phospho-Ser392) antibody (orb559175)

Description	P53 (phospho-Ser392) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p53 around the phosphorylation site of Ser392 (PD(p-S)D)
Target	phospho-P53 (Ser392)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	43 kDa
Uniprot ID	P04637
Dilution Range	WB=1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Urinary bladder (Mou...)