

Product Datasheet

Human FGF-2 protein (orb55918)

Description	Recombinant human FGF-2 protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,200mM NaCl, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Application notes	(Q65I,C96S,N111G)
Protein Sequence	Met1-Ser155
Source	E.coli
Uniprot ID	BAG70264.1
Expiration Date	6 months from date of receipt.