


Product Datasheet

PLC Gamma 1 antibody (orb559181)

Description	PLC Gamma 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Phospholipase C gamma 1 protein, around C-terminal 150aa
Target	PLCG1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	148 kDa
Uniprot ID	P19174
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-PLCG1 a...