

Product Datasheet

IRF3 (phospho-Ser386) antibody (orb559184)

Description	IRF3 (phospho-Ser386) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRF3 around the phosphorylation site of Ser386 (S(p-S)LE)
Target	phospho-IRF3 (Ser386)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	47 kDa
Uniprot ID	Q14653
Dilution Range	ICC=1:50, WB=1:500-1000, IHC-P=1:100-500, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.