



Product Datasheet

Akt (phospho-Thr450) antibody (orb559195)


Description	Akt (phospho-Thr450) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Akt around the phosphorylation site of Thr450 (IT(p-T)PP)
Target	phospho-Akt (Thr450)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	56 kDa
Uniprot ID	P31749
Dilution Range	IHC-P=1:50-100, WB=1:500-2000, IHC-F=1:50-100
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse brain
using phos...



IHC-P
analysis of
rat brain
using
phosph...



IHC-P
analysis of
Human
kidney using
pho...