


Product Datasheet

AKT (phospho-Ser124) antibody (orb559196)

Description	AKT (phospho-Ser124) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Akt1 around the phosphorylation site of Ser124 (SG(p-S)PS)
Target	phospho-AKT (Ser124)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	56 kDa
Uniprot ID	P31749
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Mouse Lung; Lane 2: Ra...