

Product Datasheet

CRISPR-Cas9 antibody (orb559217)

Description	CRISPR-Cas9 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fragment corresponding to Streptococcus pyogenes CRISPR-Cas9 (N terminal)
Target	CRISPR-Cas9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	160 kDa
Uniprot ID	Q99ZW2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.