


Product Datasheet

Bcl-XL antibody (orb559229)

Description	Bcl-XL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Bcl-XL protein, around 1-100aa
Target	Bcl-XL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	26 kDa
Uniprot ID	Q07817
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot:HeLa(Human)
Cell;Jurkat(Hum...