

Product Datasheet

Human CD10 protein (orb55932)

Description	Recombinant human CD10 protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	(N-8xHis)
Protein Sequence	Tyr52-Trp750
Source	Human Cells
Uniprot ID	P08473
Expiration Date	6 months from date of receipt.