

Product Datasheet

Human TNFRSF10D protein (orb55935)

Description	Recombinant human TNFRSF10D protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	(C-Fc)
Protein Sequence	Ala56-His211
Source	Human Cells
Uniprot ID	Q9UBN6
Expiration Date	6 months from date of receipt.