

Product Datasheet

Human Vaspin protein (orb55945)

Description	Recombinant human Vaspin protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	(C-10His)
Protein Sequence	Leu21-Lys414
Source	Human Cells
Uniprot ID	Q8IW75
Expiration Date	6 months from date of receipt.