

Product Datasheet

CNTF Receptor Alpha antibody (orb559489)

Description	CNTF Receptor Alpha antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CNTF Receptor alpha (71-170/372 aa)
Target	CNTF Receptor alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P26992
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.