

Product Datasheet

Human IL-1 R3 protein [Out of stock] (orb55952)

Description	Recombinant human IL-1 R3 protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Tag	C-Fc&His
Note	For research use only
Application notes	(C-Fc-6His)
Protein Sequence	Ser21-Gln356
Uniprot ID	Q9NPH3
Expiration Date	6 months from date of receipt.