

Product Datasheet

Mouse mGK-6 protein (orb55955)

Description	Recombinant mouse mGK-6 protein.
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	(C-6His)
Protein Sequence	Pro19-Asp261
Source	Human Cells
Uniprot ID	P15947
Expiration Date	6 months from date of receipt.