

Product Datasheet

BRLF1 antibody (orb559588)

Description	BRLF1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from HHV4tp2 BRLF1 (81-180/605 aa)
Target	BRLF1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66.5 kDa
Uniprot ID	Q3KSS7
Expiration Date	12 months from date of receipt.