

Product Datasheet

Dnmt3a (phospho-Ser390 + Ser393) antibody (PE)
[Out of stock] (orb559995)

Description

Rabbit polyclonal antivosy

+ Ser393) conjugated to PE

Species/Host

Rabbit

Conjugation

PE

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Dnmt3a around the phosphorylation site of Ser390 + Ser393

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearence

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

102kDa

Uniprot ID
[Q1LZ53](#), [Q9Y6K1](#)
Expiration Date

12 months from date of receipt.