

Product Datasheet

MAPKAPK2 (phospho-Thr334) antibody (orb560142)

Description	MAPKAPK2 (phospho-Thr334) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAPKAPK2 around the phosphorylation site of Thr334 (PQ(p-T)PL)
Target	Phospho-MAPKAPK2 (Thr334)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	P49137
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.