

Product Datasheet

CD1A antibody (Biotin) (orb560520)

Description	CD1A antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CD1A
Target	CD1A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P06126
Expiration Date	12 months from date of receipt.