

Product Datasheet

Beta Catenin antibody (Biotin) (orb560526)

Description	Beta Catenin antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombinant human β -Catenin
Target	Beta catenin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P35222
Expiration Date	12 months from date of receipt.