

Product Datasheet

Transferrin antibody (HRP) (orb561005)

Description	Transferrin antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Full length native Transferrin protein purified from human plasma
Target	Transferrin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
CAS Number	11096-37-0
Uniprot ID	P02787
Dilution Range	WB=1:500-1000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.