

Product Datasheet

**Nucleoplasmin 2 antibody (HRP) [Out of stock]
(orb561010)**

Description

Rabbit polyclonal antivosy against

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from mouse Nucleoplasmin 2

Preservatives

Constituents: 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin, 50% glycerol. Or Lyophilized. Buffer = 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin. Reconstitute with sterile distilled water.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

24kDa

Uniprot ID
Q86SE8, Q80W85
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200

Expiration Date

12 months from date of receipt.