

Product Datasheet

ATPase Na⁺ antibody (HRP) (orb561012)

Description	ATPase Na ⁺ antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ATPase Na ⁺ / K ⁺ beta 2
Target	ATPase Na ⁺ / K ⁺ beta 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14415
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.