

Product Datasheet

CD41 antibody (HRP) (orb561030)

Description	CD41 antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human ITGA2B
Target	CD41/ITGA2B
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08514
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.