

Product Datasheet

Hsp90 antibody (HRP) (orb561032)

Description	Hsp90 antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Hsp90
Target	Hsp90
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P07900
Expiration Date	12 months from date of receipt.