


Product Datasheet

Bovine BOLA-DQA1 ELISA Kit (orb561301)

Description	Bovine BOLA-DQA1(Bos taurus major histocompatibility complex class II DQ alpha 1) ELISA Kit
Reactivity	Bovine
Range	1.563-100ng/ml
Target	BOLA-DQA1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A6H6Y1
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Bovine
BOLA-
DQA1(BOLA-
DQA1) ELISA
Kit