


Product Datasheet

Bovine IGF-1 ELISA Kit (orb561321)

Description	Bovine IGF-1 ELISA Kit
Reactivity	Bovine
Range	12.5-800ng/ml
Target	IGF-1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P07455
Sensitivity	7.5ng/ml
Expiration Date	Please enquire.



Bovine IGF-1 (Insulin-like growth factor ...)