


Product Datasheet

Bovine alpha 2 Macroglobulin ELISA Kit (orb561339)

Description	Bovine alpha 2 Macroglobulin ELISA Kit
Reactivity	Bovine
Range	0.156-10µg/ml
Target	A2M
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7SIH1
Sensitivity	0.094µg/ml
Expiration Date	Please enquire.



Bovine
A2M(alpha-2-
Macroglobulin)
ELISA ...