


## Product Datasheet

### Bovine IL-17A ELISA Kit (orb561385)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Bovine IL-17A ELISA Kit |
| <b>Reactivity</b>      | Bovine                  |
| <b>Range</b>           | 15.625-1000pg/ml        |
| <b>Target</b>          | IL-17A                  |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>Q687Y7</b>           |
| <b>Sensitivity</b>     | 9.375pg/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Bovine IL-17A  
(Interleukin-17) ELISA  
Kit