


## Product Datasheet

### Gallus S100A6 ELISA Kit (orb561484)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Gallus S100A6 ELISA Kit |
| <b>Reactivity</b>      | Gallus                  |
| <b>Range</b>           | 31.25-2000pg/ml         |
| <b>Target</b>          | S100A6                  |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>Q98953</b>           |
| <b>Sensitivity</b>     | 18.75pg/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Chicken  
S100A6(S100  
Calcium  
Binding  
Prot...