


Product Datasheet

Gallus AP13 ELISA Kit (orb561497)

Description	Gallus AP13 ELISA Kit
Reactivity	Gallus
Range	31.25-2000pg/ml
Target	AP13
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q5F3U0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Chicken
AP13(Apelin
13) ELISA
Kit