


## Product Datasheet

### Gallus DPD ELISA Kit (orb561499)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Gallus DPD ELISA Kit                   |
| <b>Reactivity</b>        | Gallus                                 |
| <b>Range</b>             | 3.12-200ng/ml                          |
| <b>Target</b>            | DPD                                    |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Sensitivity</b>       | 1.875ng/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Chicken  
DPD(Deoxypyridinoline)  
ELISA Kit