


## Product Datasheet

### Guinea pig EGF ELISA Kit (orb561570)

|                        |                          |
|------------------------|--------------------------|
| <b>Description</b>     | Guinea pig EGF ELISA Kit |
| <b>Reactivity</b>      | Guinea pig               |
| <b>Range</b>           | 15.625-1000pg/ml         |
| <b>Target</b>          | EGF                      |
| <b>Storage</b>         | 2-8 °C for 6 months      |
| <b>Note</b>            | For research use only    |
| <b>Uniprot ID</b>      | <b>H0UUP9</b>            |
| <b>Sensitivity</b>     | 9.375pg/ml               |
| <b>Expiration Date</b> | Please enquire.          |



Guinea pig  
EGF(Epidermal  
Growth  
Factor) ...