


Product Datasheet

Guinea pig ET-1 ELISA Kit (orb561571)

Description	Guinea pig ET-1 ELISA Kit
Reactivity	Guinea pig
Range	1.25-80pg/ml
Target	ET-1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P97740
Sensitivity	0.625pg/ml
Expiration Date	Please enquire.



Guinea pig
ET-
1(Endothelin
1) ELISA Kit