


Product Datasheet

Guinea pig IFN-gamma ELISA Kit (orb561575)

Description	Guinea pig IFN-gamma(Interferon Gamma) ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/ml
Target	IFN- γ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8CGS0
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Guinea pig
IFN- γ (Interferon
Gamma)
ELIS...