


## Product Datasheet

### Guinea pig IL-10 ELISA Kit (orb561581)

|                        |                            |
|------------------------|----------------------------|
| <b>Description</b>     | Guinea pig IL-10 ELISA Kit |
| <b>Reactivity</b>      | Guinea pig                 |
| <b>Range</b>           | 62.5-4000pg/ml             |
| <b>Target</b>          | IL-10                      |
| <b>Storage</b>         | 2-8 °C for 6 months        |
| <b>Note</b>            | For research use only      |
| <b>Uniprot ID</b>      | <b>Q9Z1Y5</b>              |
| <b>Sensitivity</b>     | 37.5pg/ml                  |
| <b>Expiration Date</b> | Please enquire.            |



Guinea pig IL-10(Interleukin-10) ELISA K...