


## Product Datasheet

### Guinea Pig TGF beta 1 ELISA Kit (orb561606)

|                        |                                 |
|------------------------|---------------------------------|
| <b>Description</b>     | Guinea Pig TGF beta 1 ELISA Kit |
| <b>Reactivity</b>      | Guinea pig                      |
| <b>Range</b>           | 31.25-2000pg/ml                 |
| <b>Target</b>          | TGF-β1                          |
| <b>Storage</b>         | 2-8 °C for 6 months             |
| <b>Note</b>            | For research use only           |
| <b>Uniprot ID</b>      | <b>Q9Z1Y6</b>                   |
| <b>Sensitivity</b>     | 18.75pg/ml                      |
| <b>Expiration Date</b> | Please enquire.                 |



Guinea Pig TGF-β1(Transforming Growth F...