


Product Datasheet

Human BDNF ELISA Kit (orb561648)

Description	Human BDNF ELISA Kit
Reactivity	Human
Range	31.2-2000pg/ml
Target	BDNF
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	The linearity of the kit was assayed by testing samples spiked with appropriate concentration of BDNF and their serial dilutions. The results were demonstrated by the percentage of calculated concentration to the expected.
Uniprot ID	P23560
Sensitivity	2pg/ml
Expiration Date	Please enquire.



Line graph illustrates about the Ag-Ab r...