


## Product Datasheet

### Human CHI3L1 ELISA Kit (orb561697)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Human CHI3L1 ELISA Kit |
| <b>Reactivity</b>      | Human                  |
| <b>Range</b>           | 62.5-4000pg/ml         |
| <b>Target</b>          | CHI3L1                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P36222</b>          |
| <b>Sensitivity</b>     | 37.5pg/ml              |
| <b>Expiration Date</b> | Please enquire.        |



Human  
CHI3L1(Chitinase-  
3-like protein 1)...