


## Product Datasheet

### Human CXCL6 ELISA Kit (orb561709)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human CXCL6 ELISA Kit |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 31.25-2000pg/ml       |
| <b>Target</b>          | CXCL6                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P80162</b>         |
| <b>Sensitivity</b>     | 18.75pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
CXCL6(C-  
X-C motif  
chemokine  
6) ELI...