


Product Datasheet

Human ACTalpha1 ELISA Kit (orb561780)

Description	Human ACTalpha1(Actin Alpha 1, Skeletal Muscle) ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	ACTα1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P68133
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
ACTα1(Actin
Alpha 1,
Skeletal
Mus...