


Product Datasheet

Human MIP-1alpha ELISA Kit (orb561824)

Description	Human MIP-1alpha(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
Reactivity	Human
Range	23.438-1500pg/ml
Target	MIP-1 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P10147
Sensitivity	14.063pg/ml
Expiration Date	Please enquire.



Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha)