

Product Datasheet

Human TNFSF18 ELISA Kit (orb561896)

Description Human TNFSF18(Tumor necrosis factor ligand superfamily member 18) ELISA Kit

Reactivity Human

Range 0.156-10ng/ml

Target TRAIL/TNFSF18


Storage 2-8 °C for 6 months

Note For research use only

Uniprot ID **Q9UNG2**

Sensitivity 0.094ng/ml

Expiration Date Please enquire.



Human
TRAIL/TNFSF18(Tumor
necrosis facto...