


## Product Datasheet

### Human IgG ELISA Kit (orb561949)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human IgG ELISA Kit   |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 1.563-100ng/ml        |
| <b>Target</b>          | IgG                   |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Sensitivity</b>     | 0.938ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
IgG(Immunoglobulin  
G) ELISA Kit