


## Product Datasheet

### Human TPO ELISA Kit (orb561954)

|                        |                                     |
|------------------------|-------------------------------------|
| <b>Description</b>     | Human TPO(Thrombopoietin) ELISA Kit |
| <b>Reactivity</b>      | Human                               |
| <b>Range</b>           | 31.25-2000pg/ml                     |
| <b>Target</b>          | TPO                                 |
| <b>Storage</b>         | 2-8 °C for 6 months                 |
| <b>Note</b>            | For research use only               |
| <b>Uniprot ID</b>      | <b>P40225</b>                       |
| <b>Sensitivity</b>     | 18.75pg/ml                          |
| <b>Expiration Date</b> | Please enquire.                     |



Human  
TPO(Thyroid  
peroxidase)  
ELISA Kit