


Product Datasheet

Human ACA-IgG ELISA Kit (orb561969)

Description	Human ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
Reactivity	Human
Range	3.125-200ng/ml
Target	ACA-IgG
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	indirect ELISA
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Human
ACA-
IgG(Anti-
Cardiolipin
Antibody ...