


## Product Datasheet

### Human CGA ELISA Kit (orb561975)

|                        |                                     |
|------------------------|-------------------------------------|
| <b>Description</b>     | Human CGA(Chromogranin A) ELISA Kit |
| <b>Reactivity</b>      | Human                               |
| <b>Range</b>           | 0.781-50ng/ml                       |
| <b>Target</b>          | CGA                                 |
| <b>Storage</b>         | 2-8 °C for 6 months                 |
| <b>Note</b>            | For research use only               |
| <b>Uniprot ID</b>      | <b>P10645</b>                       |
| <b>Sensitivity</b>     | 0.469ng/ml                          |
| <b>Expiration Date</b> | Please enquire.                     |



Human  
CGA(Chromogranin  
A/Pancreastatin) ...